

INSTRUCTIONS FOR NITRATE + NITRITE AS N SAMPLE COLLECTION

(Please read all instructions before collecting.)

1. The Nitrate + Nitrite as N sample kit consists of:
 - Analysis request form
 - One 125mL plastic bottle containing a disk with acid preservative.
2. The following instructions describe the method of collection:
 - ❖ **Do not rinse the sample container or discard the acid preservation disk!**
The filter disk is saturated with 50% sulfuric acid (*see warning below*) and is used to preserve the sample for the analysis of Nitrate + Nitrite as N.
 - a. Flush the sample tap for five (5) minutes prior to collecting the sample.
 - b. Reduce the flow of the water and fill the sample container to the shoulder of the bottle. Replace and tighten the cap and make sure it is not leaking.
3. Please complete the Inorganic Chemical Examination of Water for Nitrate* request form.
 - *Only original bar-coded forms will be accepted. NO COPIES.***
4. Place the sample bottle and analysis request form in the mailing container and return to the ISDH as soon as possible ***within seventy-two (72) hours.***

INDIANA STATE DEPARTMENT OF HEALTH
Chemistry Laboratories
550 WEST 16th STREET, SUITE B
INDIANAPOLIS, IN 46202

5. The results will be sent within ***twenty-one (21) days*** after receiving the sample.

WARNING

The Indiana State Department of Health assumes no liability for the use of this kit.

Precautions in Sample Collection for Nitrate + Nitrite as N Analysis:

CAUTION: Do not rinse out the bottle.

CAUTION: Do not fill the bottle past the shoulder of the sample container.

CAUTION: Do not contaminate the contents of the filled bottle.

CAUTION: The filter disk in the sample bottle is saturated with 50% sulfuric acid which is very corrosive and must remain in the bottle at all times. Sulfuric acid can burn the skin. If acid gets on the skin, immediately flush the area with cold water for 3-5 minutes. All burns of any significance should be treated by a physician.