

Facility Name _____ **Ammonia Nitrogen (NH₃-N)** SM 4500 NH₃ D.

Chlorine Present Yes / No Removed by _____

Analyst _____ Analysis Date _____ Analysis Start Time _____ Analysis Stop Time _____

Calibration Standards Prepared _____ Stock Standard Lot #: _____

Check Standard Prepared _____ Stock Standard Lot #: _____

Standards and Samples Analyzed at the Same Temperature Yes / No Check Std. _____ % Rec.

Slope Check: Initial mV⁰ _____ Final mV⁰ _____ Δ mV⁰ _____ Acceptable = 58 ± 3 mV⁰ @ 20° C
or 57 ± 3 mV⁰ @ 25° C

Calibration Data

Standard	Initial Reading	Adjusted Reading	Slope

Analysis Data – Samples not analyzed on Sample Date preserved to pH < 2 with H₂SO₄ and Refrigerated

Sample ID	Sample Date	Sample Comp or Grab	Date Preserved	Preserved by	Sample Volume (mL)	Final Volume (mL)	Initial Result (mg/L)	*Dilution Factor	Report Value (mg/L)
Check Std.									
Dup.									
Spike									-----
Blank									

*Dilution Factor = $\frac{\text{Final Volume (mL)}}{\text{Sample Volume (mL)}}$

Initial Sample Result (mg/L) x Dilution Factor = Report Value (mg/L)

Spike Added (mg/L) = $\frac{\text{Concentration of Spike Added (mg/L)} \times \text{Volume of Spike Added (mL)}}{\text{Final Volume (mL) of Spiked Sample used during Analysis}}$

% Spike Recovery = $\frac{[(\text{Spiked Sample Initial Result})(\text{Final Volume}) - (\text{Sample Initial Result})(\text{Final Volume})]}{\text{Spike Added}} \times 100 = \boxed{}\%$