Facility Name:		Spectrophotometric Analysis				
Analyte	Method Used	Wavelength	nm			
Analysis Date	Analyst	Check Std. Recovery	%			

	Sample/ Standard	Sample Grab or	Date	Preserved	Sample Volume	Final Volume	Initial Result	*Dilution	Report Value
Sample ID	Date	Comp	Preserved	by	(mL)	(mL)	(mg/L)	Factor	(mg/L)
Calibration Blank								1	
Std. 1 mg/L								1	
Std. 2 mg/L								1	
Std. 3 mg/L								1	
Std. 4 mg/L								1	
Std. 5 mg/L								1	
Chk Std mg/L									
Reagent Blank								1	
Dup.									
Dup.									
Spike									

^{*}Dilution Factor = <u>Final Volume (mL)</u> Sample Volume (mL)

 $Spike \ Added \ (mg/L) = \underline{Concentration \ of \ Spike \ Added \ (mg/L) \ x \ Volume \ of \ Spike \ Added \ (mL)} \\ Final \ Volume \ (mL) \ of \ Spiked \ Sample \ used \ for \ Colorimetric \ Determination$

% Spike Recovery = [Spiked Sample Initial Concentration (mg/L) – Sample Initial Concentration (mg/L)] x 100 = _____%
Spike Added (mg/L)