

**PERTUSSIS DFA AND CULTURE**

State Form 9356 (R2/7-02)

CLIA Certified Laboratory #15D0662599

INDIANA STATE DEPARTMENT OF HEALTH

Microbiology/Rabies Laboratory

P.O. Box 7203

635 N. Barnhill Drive

Indianapolis, IN 46207-7203

Phone (317) 233-8040

Fax (317) 233-8063

PLEASE TYPE OR PRINT LEGIBLY

REQUIRED PATIENT INFORMATION				PLEASE USE BLACK INK Indiana State Department of Health is to mail report to:	
Name (Last)	(First)	Age	Sex	Facility Name	
Address				Address	
City		IN		Zip Code	City
IN		Zip Code		IN	
Attending Physician				Contact Person	
Address				Phone Number ()	Fax Number ()

SPECIMEN INFORMATION					
Date Collected:			Nasopharyngeal (Preferred): <input type="checkbox"/> Yes Other (please specify):		
Was the tube incubated prior to shipment? <input type="checkbox"/> Yes hours <input type="checkbox"/> No			Comments:		
COMUNICABLE DISEASE STATISTICS					
Date Of Onset:	Antimicrobial Therapy Begun: <input type="checkbox"/> Yes <input type="checkbox"/> No	If Yes, indicate drug:		Date Started:	Date Ended:
Pertussis Vaccination: <input type="checkbox"/> Yes <input type="checkbox"/> No		Dates:			

DO NOT WRITE BELOW THESE LINES		
LABORATORY REPORT		
PRELIMINARY REPORT:	FINAL REPORT:	
Date:		
<input type="checkbox"/> DFA test NEGATIVE	Test Results for Bordetella: <input type="checkbox"/> DFA <input type="checkbox"/> Culture	
<input type="checkbox"/> <i>Bordetella pertussis</i> <input type="checkbox"/> <i>Bordetella parapertussis</i>	<input type="checkbox"/> NEGATIVE <input type="checkbox"/> POSITIVE <input type="checkbox"/> NO GROWTH	
	<input type="checkbox"/> <i>Bordetella pertussis</i> <input type="checkbox"/> <i>Bordetella parapertussis</i>	
<input type="checkbox"/> DFA test POSITIVE		
<input type="checkbox"/> <i>Bordetella pertussis</i> <input type="checkbox"/> <i>Bordetella parapertussis</i>	<input type="checkbox"/> Unsatisfactory	
<input type="checkbox"/> Further studies in progress; final report will follow	<input type="checkbox"/> Comments	
	Lab Number:	Date Received:
	Date of Final Report:	
	<input type="checkbox"/> Copy to EPI Resource Center	

INSTRUCTIONS

Fill out completely the upper half of the request form on the reverse side. TYPE OR PRINT CLEARLY. The report will be a photocopy returned in a window envelope to the "Name and Address for Report" the patient's physician or other appropriate medical official.

SPECIMEN COLLECTION AND TRANSPORT:

1. To take the specimen, immobilize the patient's head and gently pass the swab through the nostril into the nasopharynx.

IT IS NECESSARY TO OBTAIN TWO SPECIMENS RATHER THAN A SINGLE ONE

2. *Once the specimens have been taken, prepare 4 smears (2 SEPARATE SLIDES) by rolling the swab to cover the entire circle.*

ROLL SWAB TO MAKE SMEAR ON SLIDE

3. Place the swabs in the tube containing the semisolid transport media. Be sure to immerse the swab completely in the media. Snip wire, and re-cap tightly.

LEAVE SWAB IMMERSED IN THE MEDIA FOR TRANSPORT

4. AIR-DRY the prepared slides. Replace dry smears in plastic slide holder (2 per holder).
5. RETAIN THE INOCULATED TRANSPORT MEDIA AT 35 DEGREES C. FOR 18 TO 48 HOURS until ready to transport.

If unable to incubate, hold refrigerated until ready to package and mail. RECORD THIS INFORMATION on the front of the form within the comments area.

When ready to ship, put the dried slides into the plastic slide holder and place it and the inoculated transport media into the metal container and secure the screw cap tightly.

Fold the completed patient submission form in half twice (bottom and top) and roll it around the metal container.

Insert that into the outer cardboard container and securely tighten the screw cap before mailing.

6. Be sure to give complete patient information including immunization records and antibiotic therapy. DO NOT USE KITS if expired, leaky, broken, or contaminated.
7. It is RECOMMENDED to send the specimen to the laboratory 24 hour-next day delivery. It will insure expedited transport, faster DFA results, and better recovery.

**Parasitology**

State Form 13053 (R4/5-03)
CLIA Certified Laboratory #15D0662599

INDIANA STATE DEPARTMENT OF HEALTH

Clinical Microbiology

P.O. Box 7203

635 North Barnhill Drive
Indianapolis, IN 46207-7203

Phone (317) 233-8045

Fax (317) 233-8063

Patient's Name (Last)	(First)	Age	Sex
Patient's Address		County	
Attending Physician (if not included below)		Address	

Required Specimen Information

(Please Print or Type)
Name and Address for Report

Date Collected _____

First Patient Specimen? ☐ Yes ☐ No

Type:

Purpose:

☐ Feces☐ Diagnostic☐ Release☐ _____☐ Carrier☐ Outbreak

(specify)

Other

Comments: _____

City IN Zip CodeContact
Person _____Phone
Number () _____Fax
Number () _____

Instructions are on the reverse side

LABORATORY REPORT

Do not write below this line

Lab No. _____ Date Received _____ Date of Final Report _____

☐ No trophozoites, cysts, or ova of intestinal parasites found☐ Comments _____☐ Unsatisfactory – Please resubmit for reasons below: _____☐ Copy to ISDH Epidemiology Resource Center

PARASITOLOGY

Submission of Specimens – ISDH Container No. 4A – Intestinal Parasites: Examination of fecal specimens for amoeba and other protozoa. Cryptosporidium sp. will be tested upon request.

IMPORTANT: Mailing container contains one request form, a plastic spoon, and two specimen bottles: one contains 10% Formalin and the other (PVA) Polyvinyl Alcohol fixative.

- A. Both solutions are poisonous if taken internally. Do not eat or drink solutions. Any spills should be absorbed with paper towels and discarded as trash. If fluid contacts the skin, wash with soap and water. If fluid contacts eyes, flush with water for several minutes.
- B. Kit must be used by the expiration date stamped on the outside of the mailing label. Kits exceeding the expiration date will not be tested. Submitter agencies should return expired kits to the ISDH for recycling.
- C. Write **NAME** and **SPECIMEN COLLECTION DATE** on both the **REQUEST FORM AND COLLECTION BOTTLES**. Without this complete information the specimen **WILL NOT BE TESTED**.
- D. Carefully proceed as follows.

INSTRUCTIONS

1. Fully complete the upper half of the request form. Please type or print clearly.
2. Feces must be added to EACH bottle.
3. Collect feces in a clean bedpan, plate, cardboard container, or paper or plastic cup. (DO NOT MIX URINE WITH FECES)
4. Using the plastic spoon, place a portion of feces the size of a quarter into EACH bottle.
5. Carefully break up and stir the feces in each bottle. Screw the caps on tightly and shake vigorously to further blend the specimen. Place into metal container and tighten screw caps.
6. DISCARD THE PLASTIC SPOON PROPERLY ----- DO NOT return the spoon to the laboratory.
7. Fold and wrap the completed request form *around the outside of the metal container* and place both into outer mailing tube. Tighten screw cap firmly. Mail **PROMPTLY** and unrefrigerated via the U. S. Postal Service (USPS) First Class Mail to the laboratory.
8. The maximum in-laboratory turn-around time for Parasitology specimens is 6 working days. The usual examination time is 1-2 working days after receipt. The submitter should allow sufficient transit time for the USPS to deliver results. Results will be sent as photocopies of the original report.

INDIANA STATE DEPARTMENT OF HEALTH
Clinical Microbiology
Phone (317) 233-8045
Fax (317) 233-8063

**Enteric Bacteriology**State Form 13057 (R7/5-03)
CLIA Certified Laboratory #15D0662599**INDIANA STATE DEPARTMENT OF HEALTH**

Clinical Microbiology

P.O. Box 7203

635 North Barnhill Drive
Indianapolis, IN 46207-7203

Phone (317) 233-8045

Fax (317) 233-8063

Patient's Name (Last)	(First)	Age	Sex
Patient's Address			County
Attending Physician (if not included below)		Address	

Required Specimen Information

(Please Print or Type)

Name and Address for Report

Date Collected _____

First Patient Specimen? ☐ Yes ☐ No**TYPE:****PURPOSE:**☐ Feces☐ Diagnostic☐ Release☐ Carrier☐ Contact☐ Food Service Worker☐ _____

(specify)

_____ IN _____
Zip Code

Contact Person: _____

Phone Number: () _____

Fax Number: () _____

Instructions are on the Reverse Side		LABORATORY REPORT		Do Not Write Below These Lines	
Lab No.		Date Received:		FINAL REPORT	
		Intl:		Date:	
				Intl:	

PRELIMINARY REPORT

Date _____ Initial _____

☐ *Salmonella* species isolated. Final identification report will follow.☐ *Escherichia coli* O157 isolated. Final report will follow.☐ _____

Comments:

☐ No *Salmonella*, *Shigella*, *Campylobacter* or *Escherichia coli* O157:H7 isolated☐ No Growth on Enteric Media☐ No *Clostridium* or *Bacillus* sp. isolated☐ *Clostridium* sp. isolated. Final report from the Reference Laboratory will follow☐ *Bacillus* sp. isolated. Final report from the Reference Laboratory will follow☐ *Campylobacter jejuni*☐ *Escherichia coli* O157:H7☐ *Salmonella* serotype _____☐ *Shigella sonnei*☐ *Shigella* _____ Type _____☐ Unsatisfactory – Please resubmit.☐ _____☐ Copy to ISDH Epidemiology Resource Center

ENTERIC BACTERIOLOGY

IMPORTANT INFORMATION

Submission of Specimens – ISDH Container No. 7A – Enteric Bacteriology: Examination of Fecal Specimens for enteric pathogens.

Specimens are routinely examined for *Salmonella*, *Shigella*, *Campylobacter*, and *E. coli* 0157:H7. Rectal swabs cannot be examined due to inadequate material. For other enteric pathogens and/or additional information, please call (317) 233-8045.

Before using kit, check the expiration date of the container. **Do not use kit past the expiration date or specimen will not be tested.** If you have expired kits, notify the facility from which they were acquired and request an exchange. Please return expired kits to ISDH for recycling.

The mailing container includes a request form, a spoon for handling the fecal specimen, and one specimen bottle, labeled *Cary Blair*. Do not use this kit if it does not contain the specimen bottle. Feces must be added to the preservative. Handle the bottle carefully to avoid spilling or splashing the contents. Any spills should be absorbed with paper towels and discarded as trash. If the fluid contacts the skin, wash with soap and water. If fluid contacts the eye, flush with water for several minutes. The fluid is not poisonous, but could be an irritant to sensitive individuals. Return any spilled kits to the facility from which they were obtained.

INSTRUCTIONS

7A ENTERIC CONTAINER

PLEASE READ AND FOLLOW CAREFULLY

1. Fill out the upper half of the request form completely including patient name and collection date. **TYPE OR PRINT CLEARLY.** The report will be returned to the health care provider whose address is in the space designated for **Name and Address for Report.**
2. Collect feces in a clean, dry container such as a plastic cup. **DO NOT MIX URINE WITH FECES.**
3. Open the inner container, remove the specimen bottle, and remove the screw cap. Use the plastic spoon to pick up **ONLY TWO HEAPING SPOONSFULL** of solid feces or **FOUR SPOONSFULL IF THE SPECIMEN IS LIQUID.** Place the feces into the open bottle, secure cap **VERY TIGHTLY** (write the patient name and collection date on the bottle or the specimen will not be tested). Replace the bottle into the metal container and secure the metal screw cap firmly.
4. **DISCARD THE PLASTIC SPOON AND REMAINING FECES.**
5. Fold and wrap the completed request form **around the outside** of the metal container and place into the outer mailing tube. Secure lid firmly. Mail or transport **PROMPTLY** and unrefrigerated to the laboratory.

INDIANA STATE DEPARTMENT OF HEALTH

Clinical Microbiology
Phone (317) 233-8045
Fax (317) 233-8063

**Enteric Reference Culture**

State Form 25601 (R2/5-03)

CLIA Certified Laboratory #15D0662599

Indiana State Department of Health
Clinical Microbiology
P.O. Box 7203
635 N. Barnhill Drive
Indianapolis, IN 46207-7203
(317) 233-8045

Patient Information:

Last First

Name

Address

County

Age Sex M or F Date of Onset

Attending Physician

Address

Name & Address for Report

Required Culture Information:

Source:

Date: Isolated: Submitted:

Suspected Organism

() SA () SH () OT

Send Report c/o

Phone () Fax ()

Remarks**DO NOT WRITE BELOW THIS LINE****LABORATORY REPORT**

ISDH Lab Number

Date Received

Final Report Date

By

Preliminary Report Date By

() *Escherichia coli* 0157 (H7 results pending)

() Presumptive

() Identification in Progress

() *Salmonella* serotype

Serology:

group

antigenic pattern

() *Salmonella*() *Shigella sonnei*() *Shigella* type() Negative for *Escherichia coli* 0157() *Escherichia coli* 0157:H7

()

() Copy to ISDH Epidemiology Resource Center

**ENTERIC
REFERENCE CULTURE**

INSTRUCTIONS

- 1.** Submit an 18-24 hour **PURE SINGLE COLONY CULTURE TRANSFER** on a nutrient or infusion agar slant using a screw cap tube.

DO NOT SUBMIT CULTURES IN BROTH OR AGAR PLATES

- 2.** Fill out the request form, front and back, as completely as possible. **TYPE OR PRINT LEGIBLY.** THE NAME OF THE PATIENT AND COLLECTION DATE IS REQUIRED ON BOTH THE TUBE AND FORM. On the front side check probable identification:

SA for Salmonella

SH for Shigella

OT for other

- 3.** Pack the culture securely in an approved container to meet the current transportation (D.O.T.) regulations for mailing infectious substances. Tighten screw caps.
- 4.** The report will be a photocopy of the front side of this form, and returned to the health care provider whose address is in the space designated "Name and Address for Report".

BE CERTAIN THAT YOUR ADDRESS IS COMPLETE SO IT IS RETURNED TO THE SPECIFIC SENDER.

A preliminary report will be forwarded only if there are significant preliminary findings or extended identification studies involved. FAX reports are not routinely sent.

**CULTURE INFORMATION FROM
SUBMITTING LABORATORY**

A. Identification _____

B. Serology Results _____

C. Rapid Method Kit Used ☐ YES ☐ NO

D. Type _____ Profile# _____

E. Gram Stain + - Oxidase + -

Other Information or Comments

**Indiana State Department Of Health
Clinical Microbiology
Phone (317) 233-8045
Fax (317) 233-8063**

**Reference Bacteriology Culture Identification**

SF 35898 (R3/10-01)

CLIA Certified Laboratory #15D0662599

INDIANA STATE DEPT. OF HEALTH LABORATORIES
MICROBIOLOGY LABORATORY
635 N BARNHILL DR. RM 2023
P.O. BOX 7203
INDIANAPOLIS, IN 46207-7203

REFERENCE BACTERIOLOGY LABORATORY
(317) 233-8040
Fax (317) 233-8063

☐ Mail/Fax Copy of the Report to Communicable Disease**RETURN THIS ORIGINAL FORM WITH THE SPECIMEN, NOT A COPY****REQUIRED PATIENT INFORMATION**

Name (Last, First, Middle)			
Address			
Age	Gender	Date of Onset	Physician
Diagnosis			

REQUIRED CULTURE INFORMATION

Isolation Source		Date Isolated	Date Submitted
EXAMINATION REQUESTED		Organism Suspected	
Identification	Confirmation		

SUBMITTER INFORMATION

Facility Name				
Address				
City	IN	Zip Code	Phone Number	Fax Number
Comments:				

DO NOT WRITE BELOW THIS LINE**FINAL REPORT**

Comments:		
Lab Number	Date Specimen Received	Date of Final Report

INSTRUCTIONS

The 10A container conforms to all State and Federal laws pertaining to the transport of etiologic agents. If other containers are used they must also conform to postal laws. The ISDH will reserve the right to refuse and/or discard any specimen(s) received in an inadequate or unsafe container. Submit only **PURE CULTURES** that are to be identified. Mixed cultures may not be accepted. It is best to make submissions on a low carbohydrate medium free of excess moisture. For **ANAEROBES** use stab cultures in a low carbohydrate medium or sealed chopped meat broth. **FASTIDIOUS ORGANISMS** can be sent as a heavy growth or either on blood agar or heart infusion slant.

DO NOT SEND CULTURES ON PETRI PLATES

Complete the top portion of this form. After packing the specimen in the inner container wrap the form around it, insert into the outer container, and affix the screw cap.

Check "**REFERENCE BACTERIOLOGY**" on the container address label to assure prompt delivery to the proper laboratory.

1. The Reference Bacteriology Laboratory services are available to Indiana medical facilities.
2. An effort on the part of the submitting laboratory to identify the isolate must have been made and those results made available to us upon request.

REFERENCE BACTERIOLOGY WORKSHEET

Date Received:

Date Set Up:

Lab Number:

DATE OF OBSERVATION / INITIAL

TEST READ DATE						TEST READ DATE					
TEST SET DATE						TEST SET DATE					

GRAM REACTION/MORPH.						LECITHINASE / LIPASE					
BLOOD REACTION						BILE GROWTH PY					
ATM. REQ. (O ₂ / CO ₂ / ANO ₂ / 5%O ₂)						CHOPPED MEAT DIGEST					
COLONY MORPHYLOGY						THIO GEL					
MOTILITY						IRON MILK					
OXIDASE											
CATALASE						CASEIN HYDROLYSIS					
OF FERM CTA						TYROSINE HYDROLYSIS / PIGMENT					
GLUCOSE						TINSDALE HALO					
XYLOSE						COAGULASE RABBIT PLASMA					
MANNITOL						FX 100 / BA 10 / NB 5 / PB 300					
LACTOSE						DNASE / THERMONUCLEASE					
SURCROSE						ONPG					
MALTOSE						PYR					
FRUCTOSE						LAP					
						STARCH HYDROLYSIS					
MacCONKEY						OF FERM CTA					
SS						ADONITOL					
CETRIMIDE						AMG					
SIMMONS CITRATE						AMYDALIN					
CHRISTENSEN'S UREA						ARABINOSE (L)					
NO ₃ REDUCTION/GAS						CELLOBIOSE					
NO ₂ REDUCTION/GAS						DULCITOL					
INDOLE						ERYTHRITOL					
TSI SLANT/BUTT						GALACTOSE					
TSI H ₂ S BUTT/PAP						GLYCEROL					
MR						INOSITOL					
VP						INULIN					
GELATIN HYDROLYSIS						MANNOSE					
LITMUS MILK						MELEZITOSE					
PIGMENT						MELIBIOSE					
GROWTH 25°C						METHY-MANNOSE					
GROWTH 35°C						RAFFINOSE					
GROWTH 42°C						RHAMNOSE					
GROWTH °C						RIBOSE					
ESCULIN HYDROLYSIS						SALICIN					
LYSINE DECARBOXYLASE						SORBITOL					
ARGININE DIHYDROLASE						SORBOSE					
ORNITHINE DECARBOXYLASE						STARCH					
NUTRIENT BR. 0% NaCl Growth						TREHALOSE					
NUTRIENT BR. 6% NaCl Growth						TURANOSE					
THAYER MARTIN GROWTH						XYLITOL					
LOEFFLER DIGEST/PIGMENT											
AMYLOSUCRASE											
XV REQUIREMENT											
PORPHYRIN											
ACETAMIDE											
SODIUM ACETATE											
PROPIONATE											
ANAEROBIC ALK. NO ₃						Preliminary ID by GLC Analysis / Date / Initial					
PHENYLALANINE/ MALONATE											
BACITRACIN / OPTOCHIN						Final Identification / Date / Initial					
BILE SOLUBILITY											
BILE ESCULIN											
6.5% SALT TOLERANCE											
METHYLENE BLUE HYDROLYSIS											
HIPPURATE HYDROLYSIS											
TELLURITE / TETRAZOLIUM											
6 µg VAN/ML GROWTH											

(+)positive (-)negative (A)acid (K)alkaline (pep)peptonization (±)weak,slight positive (d)delayed (S)sensitive growth
 (R)resistant growth (IR)indicator reduced (W)weak (NG)no growth (NR)no reaction (mm)millimeter



Mycobacteriology Test Request

SF 13701 (R6/12-03)
CLIA # 15D0662599

PATIENT NAME	<u>LAST</u>	<u>FIRST</u>	<u>BIRTHDATE</u>	<u>SEX</u> M F
ADDRESS				COUNTY
PHYSICIAN	MD ADDRESS			

<u>DATE COLLECTED</u>	() SPUTUM () CULTURE FOR ID	<u>SPECIMEN SOURCE</u>	<u>ANTI-MYCOBACTERIAL THERAPY</u> () NONE () STR () ISO () RIF () ETH () PZA
-----------------------	-------------------------------------	------------------------	--

<u>Sender Comments</u>	<u>Name and Address for Report</u>
------------------------	------------------------------------

DO NOT WRITE BELOW THIS LINE

<u>MICROSCOPIC EXAMINATION</u> Auramine-O Stain (400X) Acid Fast Bacteria () Not Found () Found () < 1 per field () 1-10 per field () > 10 per field () SPECIMEN UNSATISFACTORY _____ Date _____ By _____	<u>IDENTIFICATION</u> Mycobacterium () tuberculosis complex () avium complex () chelonae () kansasii () xenopi () fortuitum () gordonae () _____ () Not isolated () Culture is contaminated, more time required () Overgrown () Culture sent to CDC Date _____ By _____	<u>SIRE BACTEC SUSCEPTIBILITY</u> <table border="1"><thead><tr><th></th><th>µg/ml</th><th>S</th><th>R</th></tr></thead><tbody><tr><td>STREPTOMYCIN</td><td>2.0</td><td></td><td></td></tr><tr><td>ISONIAZID</td><td>0.1</td><td></td><td></td></tr><tr><td>RIFAMPIN</td><td>2.0</td><td></td><td></td></tr><tr><td>ETHAMBUTOL</td><td>2.5</td><td></td><td></td></tr><tr><td>PYRAZINAMIDE</td><td>100</td><td></td><td></td></tr></tbody></table> <p>S=Susceptible R=Resistant Date _____ By _____</p> <p><u>ISDH COMMENTS</u></p>		µg/ml	S	R	STREPTOMYCIN	2.0			ISONIAZID	0.1			RIFAMPIN	2.0			ETHAMBUTOL	2.5			PYRAZINAMIDE	100		
	µg/ml	S	R																							
STREPTOMYCIN	2.0																									
ISONIAZID	0.1																									
RIFAMPIN	2.0																									
ETHAMBUTOL	2.5																									
PYRAZINAMIDE	100																									
Lab No. _____ Date Received _____																										
() Copy to TB Control <input type="checkbox"/> FINAL REPORT																										

MYCOBACTERIOLOGY
(Examination for Tuberculosis)

Complete the top part of the request form. TYPE or PRINT legibly with black ink. The final report will be a photocopy of the front side only, returned to the submitting facility.

Submission of Clinical Specimens
(ISDH Container No. 6A – Tuberculosis)

1. Collect sputum early in the morning BEFORE the patient eats or drinks. It should be raised from the lungs, not saliva, and deposited directly into the furnished plastic container. Do not fill the plastic container more than half full. Wipe off any sputum from the outside of the plastic container before shipment.
2. After collection, tighten the cap to avoid breakage and leakage and place in the metal container.
3. Enclose the completed request form in the outer cardboard mailer along with the aluminum container and forward to the laboratory promptly.

Note: If gastric washings are to be submitted, you must neutralize the washings to about pH 7 within 30 minutes after collection or the organisms will die.

Submission of Cultures

1. Submit a pure culture on a tubed slant of mycobacterial culture medium, preferably Lowenstein – Jensen.
2. Pack the culture to prevent breakage and to conform to USPS and DOT regulations for “Interstate Shipment of Etiological Agents,” 49 CFR 173.196. Wrap the culture in absorbent material, place it in an inner container along with the request form and enclose securely in an outer shipping container. Do NOT send or deliver cultures grown on Petri plates; they will not be accepted.
3. The outer shipping container must be marked with the international biohazard symbol. Secure the outer container and write “MYCOBACTERIOLOGY” on the address label to assure prompt delivery to the laboratory.

Test Service Information

1. All samples received with completed paperwork are tested. If there is a problem with your submission, you will be promptly notified by telephone, mail, or both.
2. Drug susceptibility testing procedures are performed on M. tuberculosis only. Drug susceptibility testing on other species may be performed after consultation with laboratory personnel.
 - a. Bactec 460 S.I.R.E testing requires approximately 1 week too complete. This procedure is automatically performed for new patients and on patients who remain culture positive after 3 months of treatment.

Indiana State Department of Health
Clinical Microbiology
635 N. Barnhill Drive
Indianapolis, IN 46202
Phone: (317) 233-8042
Fax: (317) 233-8043



Mycology Test Request

State Form 6983 (R5/12-03)
CLIA #15D0662599

Date Received _____ ISDH Lab Number _____

Indiana State Department of Health Clinical Microbiology P.O. Box 7203 635 N. Barnhill Drive, Room 2023 Indianapolis, IN 46207-7203 (317) 233-8044 <u>SEND REPORT TO:</u>	Patient Last Name _____ First Name _____ M _____ Age _____ Sex _____
	Specimen Information () Culture for Identification Date Submitted _____ Source _____ Media _____ Suspected Organism _____ () Specimen/Sample Date Collected _____ Date Submitted _____ Type _____
Contact Person: _____ Phone Number: () _____ Ext. _____	
INITIAL REPORT Date _____ By _____ () Gram Stain Negative () Gram Stain Positive for Fungal Forms () Specimen cultured; allow 2-4 weeks or longer () Culture received; identification in progress	FINAL REPORT Date _____ By _____ () Fungus not Isolated () Specimen/Culture Unsatisfactory () Organism(s) Isolated/Identified: () Identified by Exoantigen Test () Identified by GenProbe
PROGRESS REPORT Date _____ By _____	FINAL IDENTIFICATION

() Copy to Chronic and Communicable Disease

MYCOLOGY

(Examination for Fungi)

INSTRUCTIONS

Fill out the request form as completely as possible. TYPE OR PRINT LEGIBLY. The report will be photocopy of the front side only, returned in a window envelope to the address you transcribe in the "Send Report To" box. The name or ID number of the PATIENT IS REQUIRED.

Submission of Cultures – PROVIDE OWN APPROVED MAILING CONTAINER

1. Submit a pure culture on an agar slant. Use a screw-capped tube only. PETRI DISH OR BROTH CULTURE SUBMISSIONS ARE UNSATISFACTORY AND WILL BE DISCARDED UPON RECEIPT.
2. Pack the culture in absorbent material to prevent breakage and leakage, place in an inner container assuring that it is tightly fastened. Fold and enclose the request form along with the culture securely in the outer mailing container, affix first class postage and mail promptly. All materials must meet the packaging requirements in 49 CFR 173.196 (USPS and DOT regulations for "Interstate Shipment of Etiological Agents").
3. The outer shipping container must be marked with the international biohazard symbol. Please write "MYCOLOGY" on the address label to avoid delays in delivery to the laboratory.