

Spectrophotometer – Total Nitrogen													
Total Nitrogen Results - Hach Methods 10242 and 10206, TNTplus 880 and TNTplus 835													
							(A)			(B)			
Analyst	Sample ID	Sample Date	Sample Time	Preservation Date/Time	Analysis Date	Analysis Time	Reagent TNTplus 880 0 to 16 mg/L TKN	Analysis Date	Analysis Time	Reagent TNTplus 835 0.23 to 13.50 mg/L NO ₃ -N (Nitrate/Nitrite)	Analysis Date	Analysis Time	Total Nitrogen
													Add A + B = Total Nitrogen
													Add A + B = Total Nitrogen
													Add A + B = Total Nitrogen
													Add A + B = Total Nitrogen
													Add A + B = Total Nitrogen
													Add A + B = Total Nitrogen
													Add A + B = Total Nitrogen
													Add A + B = Total Nitrogen
													Add A + B = Total Nitrogen
													Add A + B = Total Nitrogen
													Add A + B = Total Nitrogen
													Add A + B = Total Nitrogen

If the result is outside the range for the TNTplus Reagent used, it must be analyzed again using a dilution or a different range of TNTplus Reagent.

All samples must be preserved by bringing the pH below 2 (add sulfuric acid) and storing at or below 6 degrees C. **If not preserved, Nitrate and Nitrite must be analyzed separately.** Samples must be brought to pH 7 and 20 degrees C before analysis.