











FIRST AID & HEALTH

	Item	Brand	Model	Grainger Item	Features
	Hand Sanitizer, Gel, 1200 mL Cartridge	ZEP	124616	43NR34	60% Alcohol, Hand sanitizer contains agents to kill bacteria and germs that may cause illness. It can be rubbed onto hands when access to soap and water is not available. Fits with Fuzion dispenser (below)
	Fuzion Automatic Dispenser	ZEP	S93301	43NT45	Automatic wall dispensers release a premeasured amount of hand soap or sanitizer. The infrared sensors that detect the presence of hands for dispensing. These dispensers mount to a wall to free up counter and floor space in restrooms or office break rooms.
			S93306	54ZP79	
	Hand Sanitizer, Gel, 1 Gallon, 4 pack	ZEP	087824	54ZT35 (no pump)	70% Alcohol. Hand sanitizer contains agents to kill bacteria and germs that may cause illness. It can be rubbed onto hands when access to soap and water is not available. Meets current CDC & FDA guidelines as of July 2020
			087825	60JJ16 (includes pump)	
	Foot Operated Floor Stand, Fits 1 Gallon Pump Bottles	UBT	565189	60HL56	Gray, Steel, Foot Operated. Universal for any one-gallon bottle with plastic pump.
	Surface Disinfectant, Pack of 6	Purell	3342-06	49Z655	Surface Disinfectant, 32 oz Container Size, Trigger Spray Bottle Container Type, Citrus Fragrance.

FIRST AID & HEALTH

	Item	Brand	Model	Grainger Item	Features
	Hard Hat Personal ID Card	VitalID	WSID-02	45PC77	Adheres to exterior of hard hat, the ID's data window provides a name details. Water-proof sleeve stores vital medical and contact info. Non-conductive, reflective material. 100% waterproof and retains all information. Pack of 25
	Universal Plastic Personal ID Card	VitalID	WSID-05	48TD33	Strong 'Emergency Response' lid graphic. Waterproof synthetic paper ID card. Weather-proof design with lid gasket keeps recorded information safe and secure against the elements. A key-ring hook is included for easy use. Pack of 10
	Hand Sanitizer	RDC	48513501592		4 oz. Hand Sanitizer w/ Aloe, kills 99.99% of most common germs that may cause illness. Waterless Squeeze Bottle, 24 PK
	Pre-Poison Ivy Barrier, Gel	IvyX	83666	483T73	Gel, 4 oz bottle. Acts as an invisible barrier to help protect against outdoor irritants. Dries quickly, non-greasy. No wait time, simply apply to exposed skin and immediately start your outdoor task. Washes off easily with soap and water.
	Poison Oak/Ivy Cleanser	IvyX	84666	483T65	4 oz bottle, Effective in assisting the removal of poisonous plant oils (Urushiol) found in Poison Ivy, Poison Oak and Poison Sumac. How to use: Apply to all exposed areas within 1 hour of coming in contact. Clean tools and equipment before reuse. Use as laundry additive to help clean contaminated clothing.







FIRST AID & HEALTH

FIRST AID ONLY










Unitized First Aid Kit

- 25-50 people served
- Metal Case, Wall Mount
- ANSI Z308.1-2015, Class A, Type III
- (1) Fabric Bandages 1"x3" 16/Box (1) First Aid Tape (1) CPR Mask, (1) Scissors (1) Triple Antibiotic Ointment Packets 12/Box (1) Antiseptic Towelettes 10/Box (1) Burn Dressing 4"x4" (1) First Aid/Burn Cream Packets 0.9g 12/Box (1) Cold Compress 4"x5" (1) Eyewash 1 oz. w/ (2) Eye Pads & (2) Strips (1) Hand Sanitizer Packets 0.9g 10/Box (1) Nitrile Exam Gloves 4/Box (1) Conforming Gauze Roll 2" 2/Box (1) Gauze Dressing Pads 3"x3" 4/Box (2) Trauma Pad 5"x9" 2/Box (1) Triangular Sling/Bandage 40"x40"x56" (1) First Aid Guide (1) Knuckle & Fingertip Fabric Adhesive Bandages 9/Box (1) 2"x4" Patch Fabric Adhesive Bandages 6/Box (1) Povidone-Iodine Antiseptic Wipes 10/Box (1) Hydrocortisone Anti-Itch Cream Packets 12/Box (1) Sting Relief Wipes 10/Box











	Item	Model	Grainger Item	Features
	First Aid Kit	90600	48ZE22	General Purpose
	½" First Aid Tape	8-001G	39N993	Tape, White, 5 ft. L, 2 roll PK
	Antibiotic Ointment	12-001G	39N914	Triple Antibiotic Ointment, 12/box, 0.02 oz. Bacitracin Zinc, Neomycin Sulfate, Polymyxin B Sulfate
	First Aid / Burn Cream	13-006G	39N910	First Aid / Burn Cream, 0.03 oz. 0.13% Benzalkonium Chloride, 0.5% Lidocaine HCL Meets ANSI Z308.1-2015
	4"x5" Instant Cold Pack	21-004G	39P018	Cold Pack / Cold Compress 1 per box Meets ANSI Z308.1-2015
	Eye Flush & Dressing Pack	7-090G	49H429	Eye Care Kit, 2 Sets in 1 box. Eye Flush, eye covering & 2 strips. Meets ANSI Z308.1-2015



FIRST AID & HEALTH

	Item	Model	Grainger Item	Features
	Hand Sanitizer	750013	48ZV46	Instant gel hand sanitizer kills 99.9% of germs. With Vitamin E and aloe. 10 pack. Active ingredient: ethyl alcohol (62%). 0.9g per packet. Meets ANSI Z308.1-2015
	Nitrile Exam Gloves	21-026G	39P031	Powder-free, disposable. Exam/medical grade. 2 pair/box. Avoid contact with blood borne pathogens. Meets ANSI Z308.1-2015
	Conforming Gauze Roll	3JMC9	3JMC9	Sterile gauze provides form-fitting basic dressing. 2"x5 yd Roll. Use to clean wounds or wrap, or cut to length to bandage them using first aid tape. Meets ANSI Z308.1-2015
	Gauze Compress	AN205-GR	32LA69	Trauma Absorbent Gauze Compress 5"x9" 1 pad per box Meets ANSI Z308.1-2015
	Knuckle & Fingertip Bandages	1-014G	39N893	Adhesive Bandages: 4 Knuckle & 5 Fingertip. 9ct Box. Latex free. Provides coverage for knuckles to cover wounds, help stop bleeding, prevent infection, and promote healing. Meets ANSI Z308.1-2015
	2"x3" Fabric Bandages	1-007G	39N871	6 ct. Beige. Latex Free, For Use With Wound Care. Meets ANSI Z308.1-2015
	1"x3" Fabric Bandages	1-008G	39N868	16 ct. Latex-Free. Treat/protect minor cuts, abrasions to help prevent infection and promote healing. Meets ANSI Z308.1-2015
	PVP Iodine Wipes	12-015G	39N899	10 ct. Single use. Non-irritating antiseptic and germicide wipe. 50% Ethyl Alcohol 2% Iodine Meets ANSI Z308.1-2015
	Anti-Itch Cream	18-012G	39N941	12 ct. Temporary relief of itching from with minor skin irritations, inflammations and rashes. Single-use. 1% Hydrocortisone Acetate Meets ANSI Z308.1-2015

FIRST AID & HEALTH

	Item	Model	Grainger Item	Features
	Burn Dressing Packet	16-004G	39N950	Non-adherent. 4"x4" Pad. Soaked in a water-based gel. Provides cooling relief to help stop pain & further tissue damage. Helps prevent contamination/infection. Meets ANSI Z308.1-2015
	Triangular Muslin Bandage	4-006G	39N977	40"x40"x56" 1 ct. Non-woven. Can be used to dress wounds, serve as a sling or splint binder for muscular or bone injuries. Includes 2 safety pins. Meets ANSI Z308.1-2015
	CPR Mask, Tape and Scissors	90638	480G01	1 ct. Universal size. Single use, CPR Mask, First Aid Tape and Scissors. Meets ANSI Z308.1-2015
	Sting Relief Wipes	19-002G	39N932	10 ct. Provides temporary relief of pain, itching from bug bites/stings. Ethyl Alcohol 50%, Lidocaine HCL 2%. Meets ANSI Z308.1-2015
	Antiseptic Towelettes	12-018G	39N904	10 ct. Antiseptic and germicide to help prevent infection in minor cuts, scrapes, minor burns. Clean skin without the use of soap and water. Benzalkonium Chloride 0.133%. Meets ANSI Z308.1-2015
	Sterile Gauze Pad	3-0010G	39N958	4ct. Use to clean and dress wounds and stop bleeding in minor injuries or after a blood draw. Individually wrapped pads are sterile unless packaging has been damaged or opened. Meets ANSI Z308.1-2015
	First Aid Guide	BK021	11A320	Use this First Aid Only® 5-page, multifold guide to basic first aid with form to record auto accident information. Size: 3" x 4".
	Scissors	FAE-6004	65DG26	Refill Scissors, 1 1/2" x 2 1/4" x 4 1/2" ANSI Standards ANSI Z308.1-2015

FIRST AID & HEALTH

	Item	Model	Grainger Item	Features
	Tourniquet	91160	60HZ98	Black 2 in x 2 in First Aid Tourniquet
	Stop the Bleed Kit	91137	55PK45	EMS/Trauma Kit. (1) Easy Tear-open Bag, (1) Windlass Tourniquet, (1) Back-Up Elastic Tourniquet, (1) Hemostatic Blood-Clotting Gauze Pad, (1) Israeli Bandage, (2) Non-Vented Chest Seals, (1) Aluminized Rescue/Shock Blanket (1) Trauma Pad Dressing, (4) Protective Nitrile Exam Gloves, (1) Bandage Shears, (1) Marker to note tourniquet time, (1) Bleeding Control Product Usage Guideline

FIRST AID & HEALTH

BLOODBORNE PATHOGEN KIT

Unitized First Aid Kit



- Volume Absorbed per Pkg - 1 Use
- OSHA regulation 29CFR1910.1030(d)(3)(i)
- (3) BZK Antiseptic Towelettes, (1) Disposable Gown with Full Sleeves, (1) Disposable Bonnet, (2) Disposable Shoe Covers, (1) Eye Shield with Ear Loop Mask, (1) Fluid Control Solidifier Pack, 21g., (1) Biohazard Scoop, (2) 24 in x 24 in Biohazard Bags, (2) 8 in x 12 in Clear Plastic Bags, (2) Twist Ties, (3) Disposable Clean-Up Towels, (2) 5 in x 8 in Germicidal (Kills Germs) Wipes, (1) CPR One-Way Valve Face Shield, (2) Exam Quality Vinyl Gloves, 1 pair, (1) Bodily Fluid Pick-Up Guide

	Item	Model	Grainger Item	Features
Bloodborne Pathogen Kit	First Aid Only	217-O	36LF21	Plastic Case 9 in x 8-3/8 in x 2-1/2 in



BURN WRAP

- Sterile, Blue, Wool
- 30 wide x 36 Length
- FDA Approved
- Comes in a convenient pouch

	Item	Model	Grainger Item	Features
Burn Wrap	Waterjel	G3630P-4.69.000	415J32	30 in x36 in, FDA Approved

FIRST AID & HEALTH

ZOLL

Automatic Defibrillator with 5-year Program Management



- **Shock Delivered Automatically**

When the AED Plus detects a shockable heart rhythm, it delivers shock on its own, with no rescuer interaction required, no need to push a button. Once the heart analysis begins and the Fully Automatic AED Plus prompts, "DON'T TOUCH PATIENT," if a shock is needed, it will be delivered automatically after a brief verbal countdown.

- **Measuring Chest Compression with Real CPR Help**



The force required to depress the chest during CPR varies with the patient's size and build. The standard measure of chest compression quality, however, is not force but depth. The Real CPR Help® technology in ZOLL's CPR-D-padz® includes hand-placement locator, accelerometer, electronics, and a processing algorithm that work together to measure vertical displacement in space as each compression occurs.

- **Simplified Placement**

ZOLL's one-piece CPR-D-padz Electrode simplifies this step by guiding placement of the red crosshairs at the center of the imaginary line connecting the patient's nipples. Once in place, the hand-locator and the two electrode pads fall naturally into optimal position for both defibrillation and CPR.

- **Five-Year Shelf Life**

All AED electrodes transmit defibrillating electricity into the patient via metal in close contact with a salt-infused gel that is positioned between the metal and the skin. Over time, however, the salt in the gel will corrode the metal and eventually compromise electrode functionality. ZOLL's novel electrode design includes a sacrificial element that prevents significant corrosion for five years, which is unmatched in the market today.

	Item	Model	Grainger Item	Features
	Fully Automatic AED	8000-004011-01	49K147	Includes: AED, Adult pads, Batteries, Carrying Case, Rescue kit with razor scissor gloves, CPR mask, Absorbent towel, 5 year Plus Trac Management
	AED Storage Cabinet	8000-0855	4EGP6	Corrosion-resistant, highly visible defibrillator graphics protect AED from dust, dirt, and moisture. Equipped with viewing windows. 17-1/2" H x 17-1/2" W x 9" D

FIRST AID & HEALTH

ZOLL

5-year PlusTrac Program Management

PlusTrac allows you to:

- Make changes to your AED program, such as removing or adding an AED.
- Record maintenance checks of AED units.
- Run compliance reports that help manage your AED program and personnel.
- Track consumable items such as batteries and electrodes, as well as the certification expiration dates of volunteer responders.



Implementation of Users & AEDs into PlusTrac:

- Customer to provide detailed list of all expected users (assigned roles & the sites responsible for)
 - **Customer Admin (District Safety Officer)** - Oversees entire AED program - all sites & all AEDs (highest level)
 - **Site Coordinator (Maintenance Supervisor)** - Oversees specific site they are assigned to & AEDs located at that site. Multiple sites can be assigned to one person. (middle level)
 - **Inspector (SPOC/Local Contact)** - Oversees only the AEDs they are assigned to inspect. Inspectors will be responsible for logging monthly AED checks and cannot see the rest of PlusTrac or make changes (lowest level)
- All sites and users will be uploaded into PlusTrac by the customer service team
- When an AED is shipped to a site, the customer service team will reach out to the Site Coordinator listed for specific AED information
 - The site coordinator will be asked to provide:
 - AED serial number
 - Date the AED was installed in the building
 - Battery & pad expiration dates
 - Location of AED in building (i.e.: breakroom)
 - Any additional users/inspectors
- ZOLL will perform WebEx training for all assigned roles at installation. Additional individual one on one demos are available for any user needing assistance with the program, or if there is a change in roles.

FIRST AID & HEALTH

Ordering:

- Order via Grainger BPO #8631
- Please list who the Site Inspector will be for the are shipping to in the Attention To on Grainger.com.

Shipping Label / Packing List

Attention _____
Site Inspector

Setting up your AED:

- Watch the YouTube Video online: [How to set up your new AED Plus](#)
 - *Set up for Fully Automatic is the same as Semi-Automatic*
- Look for Green Check Mark on the left side of A (pictured right)

Setting up PlusTrac Management:

- Review steps on previous page for information needed
- Email tdonovan@en-proinc.com for setup instructions and support
- Attend scheduled WebEx for help in setting up

Support:

- Trainings available:
 - Reach out to your District Safety Officer for dates/scheduling
 - WebEx: How to Use AED, How to set up and use PlusTrac
 - OnSite: How to Use AED (will be conducted at DHQ)
 - [PlusTrac AED Program Management](#)
 - [PlusTrac Program Features Brochure](#)
 - YouTube Videos:
 - [Fully Automatic AED Plus Demo Video](#)
 - [AED Save Video: Saved by ZOLL AED Plus](#)
 - [ZOLL AED: Quick Overview of Features](#)
 - [Inside the Rescue: Demonstration and Practice](#)
- Need help?
 - **PlusTrac Support:** Taylor Donovan | TDonovan@en-proinc.com | 318-841-2896
 - **ZOLL South Texas Manager:** David Okeson | DOkeson@zoll.com | 281-414-6527
 - **ZOLL North Texas Manager:** Mikhail Castillo | MCastillo@zoll.com | 214-883-2019
 - **ZOLL Tech Support:** Techsupport@zoll.com | 978-421-9655
 - **Grainger AM- South, West : Christina Jones** | Christina.Jones@grainger.com | 512-202-4190
 - **Grainger AM – East Region:** Reuben Khan | Reuben.Khan@grainger.com | 713-825-6846
 - **Grainger AM- North Region:** Joshua Campbell | Joshua.Campbell@grainger.com | 214-882-0391



Any Issues:

- Any issues with your AED or a **X** appears, reach out to Zoll Tech Support immediately.
 - A replacement unit will be shipped out via Next Day Air to replace
 - A call tag will be sent to ship back original unit
 - **ZOLL Tech Support:** Techsupport@zoll.com | 978-421-9655

In the event than an AED is used:

- Call, Email, Text - support@plustrac.com | 866-352-5433
- You will receive a loner unit within 24 hours and ZOLL takes care of everything