INDIANA STATE DEPARTMENT OF HEALTH Chemistry Laboratory

Sample Collection Procedures

ARSENIC ANALYSIS IN DRINKING WATER

- 1. Obtain sample containers and submission forms from the Chemistry Laboratory. Sampling bottles must be plastic, specially cleaned for metals testing, and hold at least 125 ml.
- 2. To collect a sample:
 - a. Flush the sample tap for five (5) minutes prior to collecting a sample.
 - b. Reduce the flow of water and fill the sample container to the shoulder of the bottle.
 - c. Replace and tighten the cap and make sure it is not leaking.
- 3. Complete the analysis request form. Include your institution and email address.
- 4. Bring or Mail the water sample(s) and analysis request form(s) to:

Indiana State Department of Health Laboratories Chemistry Laboratory 550 West 16th Street, Suite B Indianapolis, IN 46202-2218

5. For Questions, please contact Mary Hagerman, mhagerma@isdh.in.gov or 317-921-5553